

Storage method	Storage temp. (°C)	F4L embalming recipe (per kg body)	F4L preservation recipe (ca. 100 liters)	Formaldehyde concentration in body (%)	Max. storage time before use	Max. use time during surgical training and dissection courses
Tank	16 to 20	80 ml F4L-A 80 ml F4L-B 5-8 ml Formaldehyde 37%	90 l water 10 l F4L-ImmA 1 l F4L-ImmB/CFC	< 0.3	No limit	Only limited by number of procedures that can be performed on one body
Freezer	-18 to -25	40 ml F4L-A 40 ml F4L-B 80 ml F4L-C	n/a	< 0.1	No limit	14 days
		No embalming (fresh-frozen body)	n/a	n/a	1 year	1 to 2 days